

permobil



## THE CITY LIFE.

With sleek good looks and high caliber performance, the C300 Corpus offers the perfect blend of stability and agility. Features include: adjustable suspension; front and rear battery access; and a seat-to-floor height of 17.5 inches with 45 degrees of tilt. No other product in its class can even compete. Small, nimble and quick, the C300 Corpus is designed for navigating busy streets and getting in and out of tight spaces. Plus, the Corpus is the world's most popular powered seating system with an ergonomic back featuring built-in, adjustable lumbar and trunk supports. The modular design of the Corpus also allows for growth and added seat functions. An ideal pairing, the C300 Corpus is ready to take you in any direction you choose. Lead the way.

**C300**  
CORPUS



## POWER FUNCTIONS



Power adjustable tilt



R-Net electronics with ICS



Power adjustable recline



Power seat height adjustment



Power adjustable legrests

## SHROUD COLORS



Crystal Blue



Vintage Wine



Metallic Green



Platinum Silver



Black Onyx



Chili Pepper Red



Popstar Pink



Titan Blue

## SPECIFICATIONS

Maximum Speed	5 MPH
Driving Range	15-19 Miles
Max User Weight	265 lbs
Ground Clearance	3"
Battery Capacity	Group 34
Width	24.25" - 28.5"
Length	43"
Recommended Battery	Group 34

Gradient Capacity	10° Up 6° Dn
Power Seat Elevator	8"
Seat Widths	17" - 19" (Adjustable)
Seat Depths	16" - 22" (Adjustable)
Backrest Heights	16" - 22"
Seat To Floor Height	17.25" - 21.5"
Electronics (standard)	R-Net
Alternative Drive Controls	Yes w/R-Net



THE POWER OF MOBILITY

**Permobil Inc.**  
**6961 Eastgate Boulevard**  
**Lebanon, Tennessee 37090**  
**800.736.0925**  
**www.permobil.com**